




ΔΙΑ•Μ  
современная лаборатория

www.dia-m.ru  
заказ on-line

## A **Complete** Mini-Vertical Electrophoresis System



### **XCell SureLock™ offers:**

- Easy set-up
- Leak-free electrophoresis
- Multiple applications

## Easy electrophoresis



Electrophoresis doesn't have to be hard.

The XCell *SureLock™* Novex Mini-Cell employs patented technology to make electrophoresis easier than ever before. With the XCell *SureLock™* Novex Mini-Cell you'll get:

- Easy set-up—there are no awkward clamps or screws to struggle with
- Leak-free electrophoresis—gels are positively locked in place
- Trouble-free runs—"incomplete runs" are eliminated, ensuring reproducible results
- One unit that meets all your electrophoresis needs—easily adapted for use with any Novex® pre-cast gel or for any blotting procedure

### Easy set-up procedure

The XCell *SureLock™* Mini-Cell can be easily set up in just three simple steps (Figure 1). The positive-locking gel tension wedge replaces clumsy clamps and devices to save

you time and effort during set-up, guaranteeing consistent electrophoresis results.


**Figure 1 – Easy set-up procedure of the XCell *SureLock™* Mini-Cell**



Drop buffer core into the lower buffer chamber of the XCell *SureLock™* Mini-Cell, insert one mini-gel in front of the buffer core and a second mini-gel or the buffer dam behind the buffer core.



Lock the gel tension wedge in place, load samples, and fill the buffer chambers with the appropriate running buffers.



Place the cell lid on the unit and you're ready to run.

## How the XCell SureLock™ Mini-Cell works

The XCell SureLock™ Mini-Cell holds gels firmly in place with a simple gel tension wedge. There is no need to struggle with awkward set-up devices such as clamps or screws. The pre-cast gels and buffer core, a U-shaped insert with platinum wire and electrode plugs, are easily placed into the mini-cell lower buffer chamber. Pushing the lever on the gel tension wedge forward into the locked position generates an even

horizontal force that seals the gel/buffer core assembly firmly into position (Figure 2). This forms two separated buffer reservoirs—the upper and lower buffer chambers. The positive locking action of the gel tension wedge ensures a trouble-free, leak-free run every time. Table 1 shows the specifications for the XCell SureLock™ Mini-Cell.

Figure 2 – Side view of the XCell SureLock™ Mini-Cell




Table 1 – Specifications of the XCell SureLock™ Mini-Cell

Specifications	
Gel Capacity:	Up to two Novex® mini-gels or one XCell II™ Blot Module
Cell Dimensions:	12.5 cm (l) x 14.4 cm (w) x 16 cm (h)
Upper Buffer Chamber Requirement:	200 ml for Novex® mini-gels
Lower Buffer Chamber Requirement:	600 ml for Novex® mini-gels
Material:	Polycarbonate
Chemical Resistance:	The XCell SureLock™ Mini-Cell is impervious to alcohol but is not compatible with chlorinated hydrocarbons (e.g. chloroform), aromatic hydrocarbons (e.g. toluene, benzene) or acetone.

## Meets all your electrophoresis needs


The XCell SureLock™ Mini-Cell can be easily converted for use in a variety of electrophoresis applications (Figure 3). The mini-cell buffer core and dam (A) allow you to electrophorese protein or nucleic acid molecules on any Novex® mini-gel\* In addition, the

XCell II™ Blot Module (B) can be inserted for blotting mini-gels right in the XCell SureLock™ Mini-Cell. There's no need to spend money purchasing individual units for each application.

Figure 3 – Multiple applications of the XCell SureLock™ Mini-Cell

### Mini-cell buffer core and dam

For analyzing protein and nucleic acid molecules on a broad range of Novex® mini-gels including PAGE, IEF, and TBE.



The buffer core and dam can be used for running any Novex® mini-gel in the 10 cm (w) x 10 cm (h) plastic cassette. The buffer core and dam are included in the XCell SureLock™ Mini-Cell (EI0001).




### XCell II™ Blot Module

For blotting mini-gels in the XCell SureLock™ Mini-Cell.



The blot module can easily be inserted in the lower buffer chamber of the mini-cell.

### NuPAGE® Novex 4-12% Bis-Tris Gel electrophoresed in the XCell SureLock™ Mini-Cell



The gel was run with MOPS SDS running buffer for 50 minutes.

Lanes 1,5,10: Bio-Rad Broad Range Standard with NuPAGE® Reducing Agent  
 Lanes 2,6,9: Mark12™ Standard  
 Lanes 3,8: *E. coli* extract with NuPAGE® Reducing Agent  
 Lanes 4,7: Turkey muscle extract with NuPAGE® Reducing Agent

### Western Blot developed with WesternBreeze® Chromogenic Kit (anti-rabbit)



Samples were run on a NuPAGE® Novex 4-12% Bis-Tris Gel (with MES Buffer) then transferred onto a nitrocellulose membrane using the XCell II™ Blot Module. NOTE: 10 pg of h-IgG was detected with clear background and no non-specific binding.


Lane 1: Multi-Mark® Standard  
 Lane 2: 10 ng of Human IgG (h-IgG)  
 Lane 3: 1 ng of h-IgG  
 Lane 4: 100 pg of h-IgG  
 Lane 5: 10 pg of h-IgG

\* See back page of this brochure for more information on Novex® mini-gels.

## Do you already have an XCell II™ Mini-Cell?

If you already have an XCell II™ Mini-Cell, you can easily upgrade it to the XCell *SureLock*™ Mini-Cell with a *SureLock* Retrofit Kit. The *SureLock* Retrofit Kit consists of a *SureLock* mini-cell lid, gel tension wedge, and molded mini-gel buffer dam (Figure 4). Simply replace the XCell II™ cell lid, front wedge, rear wedge, and square plastic dam with the components of the *SureLock* Retrofit Kit to make using the easiest mini-cell even easier.

Figure 4 – Upgrade Your XCell II™ Mini-Cell with the *SureLock* Retrofit Kit



1. Gel Tension Wedge
2. Cell Lid
3. Mini-Gel Buffer Dam

## Pick the kit that fits

The XCell *SureLock*™ Mini-Cell is available separately or with the XCell II™ Blot Module to fit your electrophoresis needs. Table 2 can help you select the product that best fits your applications.

Table 2 – XCell *SureLock*™ configurations

Product	Contents	Application
XCell <i>SureLock</i> ™ Mini-Cell	Cell for mini-gel (including buffer core, lower buffer chamber, lid, gel tension wedge, and dam)	Running Novex® mini-gels
XCell <i>SureLock</i> ™ Mini-Cell w/ XCell II™ Blot Module Kit	Cell for mini-gel and XCell II™ Blot Module	Running and blotting mini-gels
<i>SureLock</i> Retrofit Kit for XCell II™	<i>SureLock</i> Mini-Cell lid, gel tension wedge, and molded buffer dam	Upgrade XCell II™ Mini-Cell to XCell <i>SureLock</i> ™ Mini-Cell

## Simplify your electrophoresis

For the easiest and most convenient way to perform electrophoresis, choose the XCell *SureLock*™ Mini-Cell. Call Invitrogen and place your order today.

Product	Quantity	Cat. no.	Price
XCell <i>SureLock</i> ™ Mini-Cell	1	EI0001	\$365
XCell II™ Blot Module	1	EI9051	\$250
XCell <i>SureLock</i> ™ Mini-Cell and XCell II™ Blot Module	1	EI0002	\$615
<i>SureLock</i> Retrofit Kit for XCell II™ Mini-Cell	1	EI0020	\$145

# Novex® Pre-Cast Gels Meet All of Your Electrophoresis Needs

Invitrogen offers a large selection of mini-gels to suit all of your protein and nucleic acid separation needs. All Novex® mini-gels can be used in the XCell SureLock™ Mini-Cell and are guaranteed to achieve optimal resolution.

## NuPAGE® high-performance gels for protein separations

NuPAGE® Novex Gels employ a patented neutral pH electrophoresis technology that offers significant advantages over traditional Laemmli-style SDS-PAGE gel. You'll get:

- 12-month shelf life
- Excellent resolution time after time
- Efficient transfers
- Fast run times (35-50 minutes)

## Novex® Tris-Glycine gels for protein analysis

Novex® pre-cast Tris-Glycine gels are the traditional Laemmli-style gels for protein separations. With the Novex® pre-cast Tris-Glycine gels, you'll get:

- Consistent gel performance—no lot-to-lot variations
- A wide variety of single percentage and gradient gels to meet all your needs
- A choice of 8 different well formats and 2 gel thickness to meet your sample volume requirements

## Novex® IEF gels for determining isoelectric point

Novex® pre-cast vertical isoelectric focusing (IEF) gels allow you to easily and conveniently separate soluble proteins based on their isoelectric point (pI). The Novex® IEF gels provide:

- Excellent resolution for native proteins based on their isoelectric point
- Shortest IEF protocol available and no pre-focusing required
- Broad (pH 3-10) and narrow (pH 3-7) range gels to meet your needs

## More selections available

Novex® Zymogram gels for protease detection and Tricine gels for protein separation, and TBE and TBE-Urea gels for nucleic acid analysis are also available.

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7(383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vrn@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Красноярск**  
+7(923) 303-0152  
krsk@dia-m.ru

**Казань**  
+7(843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemerovo@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

