



Product Specification

Item Number	FS-A20
Product Name	PCR Film Sealer (Smart)
Dimensions(L*W*H)	36.2*17.8*29.1cm
Size/ctn	43*24*37cm
Weight	10.35kg
Compatibility	96-well/384-well PCR reaction plate, ELISA plate.
Product Features	<ul style="list-style-type: none"> • Touch screen display for easy and intuitive operation. • The sealing film is tight, no evaporation, and good consistency. • Fast and efficient, 20 seconds to complete a flat sealing film. • The parameters are flexible and adjustable, the number of times of rolling (1-9999) and the value of rolling pressure (1-3) can be programmed. • Noise-free operation, creating a comfortable working environment. • The adapter is easy to replace and can be customized according to the experimental requirements.
Storage Condition	Room Temperature
Brand	Servicebio



3188

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

